



Panorama City Medical Center Competency Validation Tool

Competency Validation Tool: **Blood Glucose Monitoring Quality Control**
Job Category: **RN, LVN, Medical Assistant**

Name _____ Unit/Dept _____ Employee # _____

Rationale for Selection: ☐ major job function ☐ high risk ☐ low volume ☐ problem prone ☐ new equipment/technology

Ask an approved validator to perform your competency validation and turn in this completed form to: _____ by: _____

Method of Validation: Direct observation Number of observations: _____	Performance (check one for each criteria)			
Medical Record # (for competency validation): _____	Able to teach others	Performed Independently	Performed with coaching	Unable to perform
1. Turn meter on / wake up the screen				
2. Enter Operator Id by scanning or using the 0-9 soft key				
3. Select and press the QC soft key from the Patient Test Screen				
4. Enter the Strip Lot Number or scan the barcode correctly				
5. Enter the QC Lot Number or scan the barcode correctly				
6. Insert the Test Strip as shown on the meter screen				
7. Mix the StatStrip Glucose Control solution gently before each use				
8. Discard first drop of control solution from the bottle to avoid contamination				
9. Perform high and low control tests				
10. Enter "NOTE" to document cleaning meter				
11. Take appropriate action for any "FAIL" QC results and able to verbalize appropriate action				
12. Enter corrective action "NOTE" for all "FAIL" QC results				
13. Dock the meter in a Data Docking Station to automatically upload stored meter data and download updated set up information				

Final Rating <input type="checkbox"/> Expert - teaches others <input type="checkbox"/> Independent – performs independently (Check One): <input type="checkbox"/> Novice* - performs with coaching <input type="checkbox"/> Not Met* - unable to perform even with coaching		* Action Plan Required
Action Plan:		Revalidation Date

Validator's signature _____ Date _____