



South Bay Kaiser Permanente Medical Center
 25825 S. Vermont Avenue
 Harbor City, CA 90710

Training and Competency Record for POCT Glucose Testing Using NOVA StatStrip

Skills Training and Evaluation	Learning Assessment/Test	Individual Information																																	
<p>At the completion of training session, the operator should be able to complete these tasks and successfully conduct POCT Glucose Testing.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 70%;"></th> <th style="width: 15%; text-align: center;">Met*</th> <th style="width: 15%; text-align: center;">Unmet</th> </tr> </thead> <tbody> <tr> <td>Describes purpose of test.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Demonstrates compliance with Standard Precautions; wears gloves.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Describes proper sample collection and application procedure.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Enter assigned Operator ID via on screen keypad or barcode scan.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Enter or scan the test strip lot number information. <i>NOTE: Keep the Stat Strip Glucose Test Strip vial tightly closed when not in use. Test strips should only be stored in the original vial at room temp.</i></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Enter the patient's ID number via on screen keypad or barcode scan.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Place the test strip into the meter and place a drop of blood on the test strip. <i>WARNING: The test strip must fill completely upon touching the blood droplet. If the test strip does not fill completely, do not touch the test strip to the blood droplet a second time. Discard the test strip and repeat the test with a new test strip.</i></td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>View test data on the screen, accept or reject the result. Remove and dispose of the test strip.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Dock the meter in a Data Docking Station to automatically upload stored meter data and download updated setup information.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Clean meter and perform QC as required.</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </tbody> </table> <p>* All skills parameters must be "met" in order to become an authorized user of this product.</p>		Met*	Unmet	Describes purpose of test.	<input type="checkbox"/>	<input type="checkbox"/>	Demonstrates compliance with Standard Precautions; wears gloves.	<input type="checkbox"/>	<input type="checkbox"/>	Describes proper sample collection and application procedure.	<input type="checkbox"/>	<input type="checkbox"/>	Enter assigned Operator ID via on screen keypad or barcode scan.	<input type="checkbox"/>	<input type="checkbox"/>	Enter or scan the test strip lot number information. <i>NOTE: Keep the Stat Strip Glucose Test Strip vial tightly closed when not in use. Test strips should only be stored in the original vial at room temp.</i>	<input type="checkbox"/>	<input type="checkbox"/>	Enter the patient's ID number via on screen keypad or barcode scan.	<input type="checkbox"/>	<input type="checkbox"/>	Place the test strip into the meter and place a drop of blood on the test strip. <i>WARNING: The test strip must fill completely upon touching the blood droplet. If the test strip does not fill completely, do not touch the test strip to the blood droplet a second time. Discard the test strip and repeat the test with a new test strip.</i>	<input type="checkbox"/>	<input type="checkbox"/>	View test data on the screen, accept or reject the result. Remove and dispose of the test strip.	<input type="checkbox"/>	<input type="checkbox"/>	Dock the meter in a Data Docking Station to automatically upload stored meter data and download updated setup information.	<input type="checkbox"/>	<input type="checkbox"/>	Clean meter and perform QC as required.	<input type="checkbox"/>	<input type="checkbox"/>	<p>1. Once opened, test strips are good for 180 days and control solutions are good for 90 days or until manufacturer's printed expiration date whichever comes first.</p> <p style="margin-left: 20px;">a) True b) False</p> <p>2. Comments (Enter Notes) are mandatory when documenting:</p> <p style="margin-left: 20px;">a) That the meter was cleaned prior to performing a QC test b) Action taken when QC results are out of range c) Action taken when critical patient results are obtained d) All of the above</p> <p>3. Control solutions and test strips must be dated when they are opened.</p> <p style="margin-left: 20px;">a) True b) False</p> <p>4. Which of the following are true about cleaning the meter?</p> <p style="margin-left: 20px;">a) Clean the meter daily before QC testing b) Record cleaning of the meter using the Comment soft key c) Use a cloth dampened with 10% bleach or 70% isopropyl alcohol or disinfectant wipe d) All of the above</p> <p>5. The test strip must fill completely upon touching the blood droplet. If the test strip does not fill completely, do not touch the test strip to the blood droplet a second time. Discard the test strip and repeat the test with a new strip.</p> <p style="margin-left: 20px;">a) True b) False</p> <p>6. You get a critical low result on your patient who is obviously not symptomatic. How do you proceed?</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>7. Mattie, an Indonesian paraplegic has a critical high blood sugar level and her husband demands to medicate with 5 units regular insulin as they do at home. How do you proceed?</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>SCORE: _____ (Passing = 80% or higher)</p>	<p>*Name: _____</p> <p>*Title: _____</p> <p>*Employee NUID#: _____</p> <p>*Unit: _____</p> <p>*Signature: _____</p> <p>*Date of Training: _____</p> <p>Training Status:</p> <p style="margin-left: 20px;"><input type="checkbox"/> New/Initial Training <input type="checkbox"/> Annual Competency</p> <p>*Required Fields</p> <p>The trainer's signature attests that the trainee has (a) successfully completed the program and scored 80% or better on the quiz, and (b) demonstrated successful skill in performing this procedure.</p> <p>_____</p> <p>Trainer's Name (PRINT)</p> <p>_____</p> <p>Trainer's Signature</p> <p>_____</p> <p>Date</p> <p>_____</p> <p>Fax completed forms to 310-517-2758 or scan document to Dennis.S.Edora@kp.org.</p>
	Met*	Unmet																																	
Describes purpose of test.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Demonstrates compliance with Standard Precautions; wears gloves.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Describes proper sample collection and application procedure.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Enter assigned Operator ID via on screen keypad or barcode scan.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Enter or scan the test strip lot number information. <i>NOTE: Keep the Stat Strip Glucose Test Strip vial tightly closed when not in use. Test strips should only be stored in the original vial at room temp.</i>	<input type="checkbox"/>	<input type="checkbox"/>																																	
Enter the patient's ID number via on screen keypad or barcode scan.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Place the test strip into the meter and place a drop of blood on the test strip. <i>WARNING: The test strip must fill completely upon touching the blood droplet. If the test strip does not fill completely, do not touch the test strip to the blood droplet a second time. Discard the test strip and repeat the test with a new test strip.</i>	<input type="checkbox"/>	<input type="checkbox"/>																																	
View test data on the screen, accept or reject the result. Remove and dispose of the test strip.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Dock the meter in a Data Docking Station to automatically upload stored meter data and download updated setup information.	<input type="checkbox"/>	<input type="checkbox"/>																																	
Clean meter and perform QC as required.	<input type="checkbox"/>	<input type="checkbox"/>																																	